

AAS Gas Purification



Salient Features Benefits:

- Gas purification panel is used to remove the vapours impurities like Moisture, CO, CO₂, Hydrocarbon from N₂O, Argon Acetylene
- Gas purification for AAs ICP-MS is optional one can also use gas distribution panel instead of gas purification panel.

Tubing / Casing Capping

- Tubing is available in different sizes like 1/2", 1/4", 1/8", 1/16"
- Tubing available in different MOC like copper, S.S., Teflon, nylon.
- Casing capping is used to cover the above tubing to maintain the esthetic look.
- Three types casing capping available, M.S. powder coated, aluminium powder coated PVC casing capping.

Additional Information:

- **Minimum Order Quantity:** 1 per piece