



## Vortex / Cyclomixers

## **Vortex Mixers**

- Designed for mixing liquids in Schools, Laboratories & Factories
- Touch/ Continuous Operation mode Selection through bi-directional Switch
- Speed Regulation through knob provided on the control panel
- Interchangeable mixing heads for use with variety of tubes

REMICM 101 rus (2)   D	
Gran room	

**CM** - 101 Plus

Technical Data		
Models	CM - 101 PLUS	
Power [W]	60	
Shaking Movement	Orbital	
Orbital Diameter [mm]	4	
Motor Type	Shaded-Pole Motor	
Motor Rating Input / Output [W]	58 / 10	
Permissible ON time	100% power 30mins	
Speed range [rpm]	0-2500	
Run Type	Continuous / touch operation	
Dimensions (mm)	127 x 130 x160	
Weight [kg]	3.5	
Protection Class acc. To DIN 60529	IP 21	

Supply 220/240 Volts, 50 Hz

Interchangeable Mixing Heads		
Models Description		
CM-01	Universal top plate 100 mm-plate attachment	
CM-01/01	Tube adapter for 18 holes test tubes ø 10mm (use with CM-01)	
CM-01/02	Platform pad for <ø99mm tubes and small vessels (use with CM-01)	
CM-01/03	Tube adapter for 8 holes test tubes ø 20mm (use with CM-01)	



CM - 01



CM - 01/03



CM - 01/01







