



# Laboratory Centrifuges



C - 852



R - 303



C - 854 / 8

## Medico / Doctor Centrifuges

These are low priced units meeting the majority of the small capacity centrifuging requirements of medical practitioners, pathological laboratories and hospitals for

routine centrifuging tests. The units are fitted with resiliently mounted motor for vibration free performance having built in 5 step speed regulator. Available with Swing out or Angle rotor.

### Technical Data

Models		C-852	C-854/4	C-854/6	C-854/8	R-303	R-304
Capacity	ml	4 X 15	4 X 15	6 X 15	8 X 15	8 X 15	24 X 5
Type of Head		Swing Out				Angle	
No. of Tubes		4	4	6	8	8	24
Max. Speed	RPM	3500	3500	3500	3500	3800	7000
Max. RCF	'g'	1750	1850	1600	1600	1450	4760
W X D X H	mm	295 x 295 x 25	310 x 310 x 295			210 x 210 x 250	

Supply: 220 - 240 Volts 50 Hz Single Phase.

## Mini Centrifuge RM - 02

Mini Centrifuge permits scientists to spin down small samples from 0.2 ml to 1.5 ml. Blinking pilot lamp (one blink per second) for accurate time counting in seconds to facilitate short spin.

This centrifuge is ideal for Applications in PCR, Cell Separation, and Biotechnology work and in clinical labs.

### Salient Features:

- PCR Strip Rotor to accommodate 2 strips of 8 x 0.2 ml
- 6 x 1.5 ml fixed Angle Rotor
- Adapter for 0.2 ml & 0.5 ml tubes
- Fix speed - 6000 RPM
- Safety Lid interlock to prevent cover opening during centrifugation
- Quick acceleration to full speed
- **Unit supplied with PCR Strip Rotor, fixed Angle Rotor and Adapters**

Supply: 220-240 Volts 50 Hz Single Phase.

