

INTRODUCING

**WORLD CLASS**

**PREMIUM RANGE**



**RS - 24 Plus**

These Shakers are ideal for mixing and development of cultures, chemicals, solvents, and assays etc. in Microbiological, Cell Culture & Life Science laboratories.

The Unit has base assembly fabricated from heavy mild steel sections. Brushless Induction drive motor with frequency drive makes the unit suitable for continuous non-stop operation. Step less electronic frequency control ensures gentle start and maintains preset speed. Compact counter balanced drive mechanism ensures high stability and reliability even in continuous operation & uneven load distribution.

## Salient Features

- 4" Large LCD display
- Universal Platform with rounded corners to accommodate interchangeable clamps of assorted sizes
- Brushless induction motor with VFD speed control
- Elegant Plastic front
- Variable speed 20 -300 RPM
- Shaking amplitude 40 mm
- Auto start after power resumption

### Technical Data

Models	Platform Size	Max. Volume	Max. No. of clamps (to be ordered separately)					
			2000 ml.	1000 ml.	500 ml.	250 ml.	150 ml.	100ml/ 50ml
RS-24 Plus	24" x 24"	18 Litres	9	16	25	36	49	64

Supply: 220-240 Volts 50 Hz Single Phase.