



Deep Freezer



These Deep Freezers are designed for quick freezing and storing of Blood Components, Serum, Vaccines, Biological & Medical Specimens, Clinical Samples, etc. at low temperatures.

Engineered for just a little bit more, the new REMI series of Chest Freezers & Upright Freezers gives you upto -25°C instead of the conventional - 20° C. Equipped with digital display and alarm features these freezers give security needed for safe sample storage even at low temperatures.



Salient Features:

- Plastic inner chamber ; Outer body pre painted G.I.
- Quick freezing function
- Sealed drawers - prevent frost loss from open door
- Solid state digital temperature controller with alarms
- Adjustable baskets / shelves
- Low energy consumption and noise levels

Technical Data

| Model No. | Chest Freezers | | | Upright Freezers | |
|--------------------------|----------------|-----------|-----------|------------------|-----------|
| | RFC - 100 | RFC - 200 | RFC - 280 | RFV - 245 | RFV - 340 |
| Capacity (Litres) | 108 | 200 | 280 | 245 | 340 |
| Temperature Range °C | +8 °C ~ -25 °C | | | -16 °C ~ -24 °C | |
| Door | 1 | 1 | 1 | 1 | 1 |
| Baskets / Shelves | 1 | 1 | 2 | 6 | 7 |
| External Dimensions (mm) | | | | | |
| Width | 610 | 864 | 1092 | 533 | 610 |
| Depth | 508 | 533 | 610 | 610 | 635 |
| Height | 838 | 838 | 838 | 1448 | 1854 |

Supply: 220-240 Volts 50 Hz Single Phase.

